

## **CHECKLIST FOR RECEIPT OF FROZEN CELLULAR PRODUCTS**

| Name Date:  |                             |            |  |
|---|-----------------------------|------------|--|
| Date of Birth: BB number  |                             |            |  |
| Hospital number Product   |                             |            |  |
| Manufacturer/Sponsor Trial name Su  Please initial when complete. Do not leave blank fields, use N  Manufacturers product specific instructions must be checked prior to receipt  |                             |            |  |
| DATE OF RECEIPT: TIME:  |                             |            |  |
| Lot Number/Batch:   | Delivered by:  Received by: |            |  |
| Dewar Serial Number:  |                             |            |  |
| Confirmatory Checks   | Performed by                | Checked by |  |
| Waybill signed and copy retained  |                             |            |  |
| Delivery details on shipper label or waybill checked  |                             |            |  |
| Tamper-evident ties intact? Yes/No Number/s:  |                             |            |  |
| Online Portal or Logger temperature/status: Alert? Yes/No<br>Temperature on Receipt: Shipper Serial Number  |                             |            |  |
| Product Checks  |                             |            |  |
| Ensure tray of dry ice prepared before removing cells. Check specific manufopening the shipper.  Remove and check each unit separately on dry ice to verify condition and p designated storage space in LN2 vapour immediately. |                             |            |  |
| Product checked for integrity Photographs taken Saved in shared area and /Printed   |                             |            |  |
| Label on product checked to match the paperwork from supplier and internal patient information  Time out of shipper:  |                             |            |  |
| Product paperwork completed: DIN:   |                             |            |  |
| Unit/Batch Numbers:   |                             |            |  |
| Expiry Date/s: Volume/s:  |                             |            |  |
| Dose/s:   |                             |            |  |



| Storage  | •                |                   |                   |                        | Input<br>by | Checked<br>by |
|--|------------------|-------------------|-------------------|------------------------|-------------|---------------|
| Time into Storage Tank:                          |                  |                   |                   |                        |             | l by          |
| Tank   | Bag No.          | Section           | Rack              | Slot                   |             |               |
| ·  | <b>Dag</b> 1101  |                   | Mach              | 0.01                   |             |               |
|  |                  |                   |                   |                        |             |               |
|  |                  |                   |                   |                        |             |               |
| Ctorogo  | oard or invent   | on log complet    | end.              |                        |             |               |
| Storage  | card or invent   | ory log complet   | .eu               |                        |             |               |
| Temperature log printed/saved as PDF and checked |                  |                   |                   |                        |             |               |
| Docum  | entation         |                   |                   |                        |             |               |
| Docume   | entation as app  | olicable receive  | d and checked     |                        |             |               |
|  |                  |                   | QP release/C      | ertificate of Analysis |             |               |
|  |                  |                   |                   | Patient Details        |             |               |
|  |                  |                   | Saralagy racu     |                        |             |               |
|  |                  |                   | Serology resu     | lts/ Donor clearance   |             |               |
|  |                  |                   |                   | Dose/s                 |             |               |
|  |                  |                   |                   | Volume/s               |             |               |
|  |                  |                   |                   | Expiry date            |             |               |
| Manufa   | cturer/ Trial On | line Portal upd   | ated (if applicab | le)                    |             |               |
| Receipt  | forms and oth    | er paperwork p    | rovided by Manı   | ufacturer completed    |             |               |
|  | rned where ap    |                   |                   |                        |             |               |
| Reques   |                  | IIE LIIVIS        |                   |                        | T           |               |
|  | details and re   | sults input       |                   |                        |             |               |
|  |                  |                   | ts input/checked  | I in TDBB              |             |               |
|  | and ATIMPs       |                   |                   |                        |             |               |
|  |                  | d in original red | cord              |                        |             |               |
| COMME  |                  |                   |                   |                        |             |               |
| 2  | <b></b>          |                   |                   |                        |             |               |
|  |                  |                   |                   |                        |             |               |
|  |                  |                   |                   |                        |             |               |